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*owner's manual*

## EQUIPOS

*delivery units*

**DABI ATLANTE**

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## Contents

Legend.....	16
Delivery unit lateral controls .....	17
Footswitch .....	17
Installation.....	17
Electronic controls.....	17
Procedures prior to equipment use/reutilization....	17
Operation.....	17
Positioning.....	17
High-speed and Low-speed handpiece terminal.....	17
Three-way syringe use .....	18
Fiber optics use.....	18
Chair activation by the delivery unit panel.....	18
Light activation by the electronic panel.....	18
Ultrasound use .....	18
Sodium bicarbonate jet.....	19
Chip blower and water lock .....	20
Flush system use .....	20
X-Ray viewer .....	20
Cuspidor water supply activation .....	20
Cup holder water supply activation .....	20
Ultraled use.....	20
Cleaning .....	21
Maintenance.....	22
Flush system, handpieces and jet/ultrasound tank supply (Cuspidor) .....	22
Recommendations and care.....	23
Procedure in case equipment functioning is altered .....	23
Sensitivity to predictable environmental conditions in normal working situations .....	23
Procedures in case of equipment disposal .....	23
Transportation and storage conditions .....	23
After installation.....	23
Preparation prior to delivery .....	23
Dabi Atlante Authorized Service Network .....	24
Symbols.....	24
Final recommendations.....	24
Classification and type .....	24
Equipment Warranty .....	24
Delivery unit technical features (for all models)...	25
Curing light technical features .....	25
Ultrasound technical features .....	25

Troubleshooting .....	25
Troubleshooting: Delivery Units.....	25
Troubleshooting: Curing Light.....	26
Troubleshooting: Sodium Bicarbonate Jet .....	26
Troubleshooting: Ultrasound.....	26

## Presentation

Dear Client,

Welcome!

We are really proud of your preference for our products. Your familiarization with the **Dabi Atlante Delivery Unit** will help maximize their performance and increase their durability and benefits. Therefore, before starting your equipment, be sure to read this manual very attentively.

Important

All the information, illustrations and specifications contained in this manual are based on data existing at the time of its publication. We reserve the right to make any modifications deemed necessary both on the product and in this manual without previous notice.

## Legend

**TechnoAir** (fig. A) and **Techno Cart** (fig. B) **Delivery Units**

- ① Instrument holder (optional)
- ② Tray
- ③ Tray holder
- ④ Three-way syringe
- ⑤ High-speed handpiece terminal
- ⑥ Micromotor terminal
- ⑦ Sodium bicarbonate jet (optional)
- ⑧ Ultrasound (optional)
- ⑨ Arm
- ⑩ X-ray viewer (optional)
- ⑪ Electronic controls (optional)
- ⑫ Handle
- ⑬ Handpiece support
- ⑭ Ultrasound led (optional)
- ⑮ Arm lock
- ⑯ Pedestal
- ⑰ Caster

**Hasteflex Delivery Unit** (fig. C)

- ① Ultrasound (optional)
- ② Sodium bicarbonate jet (optional)
- ③ High-speed handpiece terminal

- ④ Micromotor terminal
- ⑤ Straight syringe
- ⑥ Electronic controls (optional)
- ⑦ Lateral controls
- ⑧ Arm lock
- ⑨ X-ray viewer (optional)
- ⑩ Handle
- ⑪ X-ray viewer activation (optional)
- ⑫ Instrument holder (optional)
- ⑬ Tray
- ⑭ Tray holder
- ⑮ Sodium bicarbonate tank (optional)

The serial delivery unit version features three terminals (three-way syringe, high- and low-speed handpieces).

It is possible to configure the Delivery Units with up to 5 terminals, complementing them with the ultrasound, sodium bicarbonate jet, a second terminal for a high-speed handpiece or fiber optics.

## Electronic controls (fig. D1 and D2)

### Techno Air, Techno Cart and Hasteflex Delivery Units

- ① Programmable memories
- ② Seat ascent
- ③ Backrest descent
- ④ Light
- ⑤ Spitting position
- ⑥ Cuspidor water
- ⑦ Cup holder water
- ⑧ Seat descent
- ⑨ Backrest ascent
- ⑩ Zero position

## Delivery unit lateral controls (fig. E and fig. F)

### Techno Air, Techno Cart and Hasteflex Delivery Units

- ① Arm lock
- ② X-ray viewer activation (optional)
- ③ Flush system activation (optional)
- ④ Fiber optic light activation (optional)
- ⑤ Curing Light master switch (optional) or Second Fiber Optic light activation (optional)
- ⑥ Ultrasound master switch (optional)
- ⑦ Electric Micromotor Speed selection key (optional)
- ⑧ Inversion of the Micromotor rotation (optional)
- ⑨ Ultrasound power selection key (optional)
- ⑩ Ultrasound irrigation solution outflow selection key (optional)

- ⑪ Sodium Bicarbonate Irrigation outflow selection key (optional)

## Footswitch (fig. G)

- ① Chip blower (optional)
- ② Water lock (optional)
- ③ Progressive activation footswitch

## Installation

This equipment must be unpacked and installed by an authorized Dabi Atlante technician; otherwise, warranty will be voided.

## Procedures prior to equipment use/reutilization

Prior to equipment use/reutilization, follow the operation and cleaning instructions in this manual.

## Operation

### Turning the equipment on/off

Turn on the dental chair's master switch. All equipment's functions will be enabled. Turn it off after use.

The master switch has an internal LED that remains on/off when the chair is on/off.

## Positioning

By keeping the "Arm Lock" key (fig. A, item 15 and fig. C item 8) pressed, place the delivery unit in the desired position, holding it by the handle. Release the key to fix it in such position.

## High-speed and Low-speed handpiece terminal


After connecting the high-speed handpiece, remove it from its support and activate it by pressing the footswitch.

Power (air supply) can be controlled by the operator by applying higher or lower pressure on the footswitch.

## High-speed and low speed handpiece spray adjustment

Adjustment is made by a meter located on the hose next to the terminal. Turn it clockwise to increase spray and counterclockwise to decrease it.

## Three-way syringe use

Press the green button for water outflow, the blue button for air or both buttons simultaneously for spray.

## Fiber optics (optional)

Connect the terminal (Midwest or Quick Connection) to the hand piece hose and turn on the fiber optics light activation key (fig. E and F, item 4 or 5). Remove the handpiece from its support and activate the footswitch. Air supply can be controlled by the operator by applying higher or lower pressure on the footswitch.

⚠ Turn off the fiber optics light activation key when the handpiece is not being used.

## Chair activation by the delivery unit panel (optional)

### Techno and Hasteflex models with an electronic panel

#### Sound indication configuration

1. Place the Chair at the zero position, by pressing the zero-position key (fig. D1 and D2, item 10);

2. Keep the zero-position key pressed. A long sound signal will be emitted, indicating the setting mode. Release the key after:

- The first tone (single) for configuration without a sound signal;
- The second tone (double) for configuration with a standard sound signal;
- The third tone (triple) for configuration with a polyphonic sound signal.

⚠ Regardless of the configuration selected, the work position controls will be followed by a sound signal, and the movements activated by the joystick will only emit a sound signal when the seat and backrest reach their highest or lowest position.

#### Seat movement

Press key  /  to raise or lower the seat.

#### Backrest movement

Press key  /  to raise or lower the backrest.

#### Automatic movement to the zero position

Press key  for the seat and backrest to simultaneously return to the zero position. This is the most comfortable position for the patient to get on/off the chair.

#### Work positions programming

1. Place the chair at the desired work position

using the joystick;

2. Keep keys 1, 2 or 3 pressed to save the chair's position in memory 1, 2 or 3, respectively.

⚠ It is possible to save the light's intensity and the work position at the same time, by setting it to the desired intensity prior to beginning the saving process.

#### Automatic movement to the work positions

Press and release key:

- **1** to move the chair to work position 1;
- **2** to move the chair to work position 2;
- **3** to move the chair to work position 3.

After releasing the key, the chair's seat and backrest will simultaneously move to the selected work position and the light will be set to the programmed intensity.

In order to interrupt the movement, activate any other control.

#### Automatic movement to the spitting position

In order to move the backrest to the spitting position, quickly press key .

In order to return to the previous position, quickly touch key  twice.

In order to interrupt the movement, quickly touch the spitting position key while the movement is being performed.

## Light activation by the electronic panel (optional)


Press and release key  to turn the light on or off.

⚠ The reflector must be connected to the chair's circuit.

In order to change light intensity keep the key pressed until the desired intensity is achieved.


## Ultrasound (optional)

Fit the protection cover onto the handpiece.

**Endo / Perio / Retro Tips**

Insert the selected tip into the wrench and then screw it onto the handpiece.


**Dentistry Tips**

Screw the adaptor onto the handpiece and tighten it firmly with the aid of an installation wrench.

Screw the nut a little loosely and insert the desired Dentistry Tip, pushing it flush with the stop. The set must be like in.

Tighten the nut with the installation wrench.

To change the Tip, loosen the nut, remove the Tip, insert the desired Tip and tighten the nut again.

**⚠** Never use the Dentistry Tips without the correct cooling.

Place the irrigating solution in the pet tank (refer to "Solutions used"). Turn on the Ultrasound master switch (fig. E and F, item 6) and select the power using the selection key (fig. E and F, item 9) as needed:

Indicator	I	II	III	IV
Power	25%	50%	75%	100%

Select the irrigating solution outflow using the selection key (fig. E and F, item 10). Next, remove the handpiece from its support and press the footswitch to activate it.

When the procedure has been completed, release the footswitch and place the handpiece back on its support.

**Recommendations**

The ultrasound *tips* must be unscrewed at least once a day in order to prevent the thread from jamming.

Whenever the ultrasound tips are used without cooling (water), the equipment should function for a maximum period of 4 minutes followed by an 8 minute pause in order to prevent irreparable damage to the ultrasound handpieces.

The use ultrasound is contraindicated for patients and dental surgeons with a cardiac pacemaker.

**Sodium bicarbonate jet (optional)**

To fill the sodium bicarbonate reservoir in Techno Air and Techno cart delivery units, raise the tray holder (fig. A and B, item 3).

In the Hasteflex Delivery Unit, the reservoir is located in the equipment superior part (fig. C, item 15).

**Sodium bicarbonate jet use**

Remove the sodium bicarbonate reservoir cover and if there are residues on the bottom, clean with the aspirator.

Rub the Dabi Atlante sodium bicarbonate packet for a while to homogenize it.

Pour one Dabi Atlante sodium bicarbonate packet inside the reservoir. Never pour more than the indicated by the volume level mark.

Replace the cover. Make sure the reservoir is tightly closed.

**⚠** Use only Dabi Atlante sodium bicarbonate, since other products may have varied grain size and uncontrolled moisture, which can clog the equipment.

Fill up the Pet Reservoir in the Cuspidor with your preferred fluid (check section "Reservoir replenishment").

Set the irrigant flow rate in the Selection Key (fig. E or F, item 11), remove the handpiece from the support and press foot control to start the procedure.

Upon finishing the procedure, release the foot control and place the handpiece into the holder.

**Recommendations**

If the sodium bicarbonate jet tip becomes clogged, remove the hose and unclog it.

Sodium bicarbonate dehydrates lips. Apply some type moisturizer on the patient's lips prior to any procedures.

The use of the sodium bicarbonate jet is contraindicated for patients with renal disorders or who are being submitted to hemodialysis.

Sodium bicarbonate spray may interfere with respiratory disorders. Medical advice is indicated in such cases.

The sodium bicarbonate jet handpiece must be unscrewed at least once a week in order to prevent the thread from jamming.

### Chip blower and water lock (optional)

To activate the chip blower (fig. G, item 1), press its key using your foot. Once it has been activated, air will come out through the high-speed handpiece without the need to activate the turbine and without activating spray. To interrupt, remove your foot from the key.

To activate the water lock (fig. G, item 2), set the key towards the chip blower to unblock water from the handpieces. Return it to the original position to block it.

### Flush system (optional)

Remove the handpieces from their supports and hold them by keeping the spray orifices towards the inside of the cuspidor's basin.

Activate the flush system key (fig. E and F, item 3) for 30 to 60 seconds.

Thus, the solution in the flush tank will go through the tubes and disinfect them.

After disinfection, activate the spray for at least 15 seconds in order to eliminate disinfectant residues from inside the handpieces tubes.

### X-Ray viewer (optional)

Turn the key on/off to turn the x-ray viewer on/off (fig. E and F, item 2).

### Cuspidor water supply activation (optional)

#### Techno and Hasteflex Delivery Units with an electronic panel

#### Automatic activation of the cuspidor water supply

Press key  with a quick touch in order to open or close water flow into the cuspidor.

Press key  with two quick touches to open the water flow into the cuspidor for the programmed time.

#### Time programming for water flow into the cuspidor's basin

Keep key  pressed and a polyphonic sound will

be emitted, indicating the program mode. Release the key after:

- The emission of the first tone to program for 10 seconds;
- The emission of the second tone to program for 20 seconds, successively following the same procedure for up to 60 seconds.

#### Continuous flow

For a continuous flow of water into the basin, press the cuspidor water key (fig. D, item 9) until two short tones are heard, and then release the key. To stop that command, press the key again.

This time will be programmed for 1 minute and 30 seconds although it can also be programmed for 0 second to 4 minutes. To reprogram, press the cuspidor water key (fig. D, item 9) until two short tones and one long tone are heard and release the key. When the desired time is achieved, press the key again.

#### Other times

To program water flow into the basin for 1 minute, press the cuspidor water key (fig.D, item 9) two consecutive times. For a two-minute flow, press the key three consecutive times, and successively follow the same procedure for up to 4 minutes of water flow. To stop that command before the time has expired, press the key again.


### Cup holder water supply activation (optional)

Keep key  pressed for water supply to the cup filler. Release it to interrupt the water supply.

### Ultraled (optional)

#### Curing

Fit the acrylic curing tip into the front orifice of the handpiece.

Turn the Ultraled master switch (fig. E and F, item 5).

 The acrylic curing tip is supplied with **Ultraled**. It is possible to use the semi-oriented tip, which is sold separately.

After fitting the tip into the equipment, clean it with some gauze, since oiliness from the operator's hands and other stains on the acrylic tip interfere with light conduction.

Turn on the equipment master switch.

Place the handpiece into the patient's mouth and position the tip at a maximum distance of 2 mm from the material to be cured.

Press the button to activate the LED, which will go on gradually. As soon as it is activated, the LED will remain on during the programmed time, and a tone will sound every 10 seconds.

When the set time has expired, the LED will be deactivated (refer to "Time setting").

**!** In case it is necessary to repeat the curing procedure, press the **Ultraled's** activation button again.

In order to interrupt the curing cycle, press the activation button one more time.

### Bleaching

Fit the whitening acrylic tip (1 or 3 teeth) into the front orifice of the handpiece.

**!** Whitening tips are sold separately.

**!** After fitting the tip into the equipment, clean it with some gauze, since oiliness from the operator's hands and other stains on the acrylic tip interfere with light conduction.

Turn on the Ultraled master switch (fig. E and F, item 5).

Follow the gel manufacturer's instruction concerning activation time and the distance between the tip and the gel.

Press the button to activate the LED.

The LED will go on gradually, and as soon as it is activated, it will remain on during the programmed time. A tone will sound every 10 seconds.

When the set time has expired, the LED will be deactivated (refer to "Time setting").

**!** In order to interrupt the whitening cycle, press the activation button one more time.

### Time setting

Turn the equipment's main switch on. Press the handpiece's activation button/trigger and keep it pressed until a long tone is heard.

Release the activation button. Press it the number of times to be programmed. Each touch corresponds to 10 seconds. The maximum time programmed is 100 seconds.

### Intensity time adjustment program

Press the handpiece's activation button to activate the LED. Next, without releasing the activation button, wait 3 seconds, and a long tone will be emitted. Still keeping the activation button pressed,

wait until the second long tone is emitted. Release the activation button and wait until you hear another tone, which indicates the first program phase - "Low Power". Next, press the activation button to program the application percentage. Each touch corresponds to 10% of the total pre-set functioning time while the maximum is 100% (in case the activation key is not pressed, the system will understand that no percentage has been programmed in this phase).

After the desired application percentage has been programmed, and if it is less than 100%, wait 2 seconds, and two tones will be emitted, indicating the second program phase - "Medium Power". Repeat the procedure used for the first phase, and only the remaining percentage will be allowed to be set. "High Power" will be the percentage remaining from the two first phases. Upon conclusion of the second phase and after 2 seconds, 3 fast tones will be emitted, which indicates the end of the program.

In case you want to program only some of the phases, enter the program module as explained above and do not press the activation button for 2 seconds in each phase. In this way, the system will automatically shift to the next phase. In such case, if the activation button is not pressed during the programming of phases 1 and 2, application will be at 100% of "High Power".

Example: For a 60-second application time and the phases respectively programmed at 10, 20 and 70%, the LED will remain on for 6 seconds at low power, for 12 seconds at medium power and for 42 seconds at high power.

### Cleaning

Clean the equipment surface with a clean soft cloth moistened in Aplic Odonto by Dabi Atlante or another product with similar chemical properties.




**!** The use of other chemical products is not recommended as it may damage the equipment.

### Handpiece support cleaning (Techno Air and Techno Cart models)

Remove the handpieces from their supports and place them on the tray.

Hold the cover and pull the locks with your finger to release it (fig. S, T and U).

Wash the cover using a soft sponge, water and neutral soap.

Place the cover on the support and inversely follow the procedure.

Replace the handpiece on the support.

### Ultrasound handpiece cleaning

To clean the handpiece surface, use a clean cloth moistened in Aplic Odonto®. The handpiece cover, *tips* and *tip* installation key must be washed in running water and sterilized by autoclave at 135°C.

### Sodium bicarbonate jet handpiece cleaning

Sterilize the sodium bicarbonate jet handpieces, *tips*, *tip* installation key and handpiece protection cover by autoclave at 135°C.

⚠ Remove the hose from the sodium bicarbonate jet handpiece tip once it is not autoclavable.

### Curing light cleaning

Clean the equipment surface with a clean soft cloth moistened in Aplic Odonto by Dabi Atlante or another product with similar chemical properties.

⚠ The use of other chemical products is not recommended as it may damage the equipment.

⚠ In case PVC film is used, make sure that the ventilation orifices (air inlet and outlet) have not been obstructed.  
Never wrap the acrylic tip with any type of material.

The acrylic curing tip and the acrylic whitening tips must not be autoclaved.

### Syringe Cleaning

Completely remove any barrier placed on the syringe (PVC film, for instance). Next, remove the syringe tip, clean it and send it for autoclave sterilization.

The syringe's body must be disinfected according to the disinfection protocol used in the dental office. Thoroughly clean the syringe's surface using detergent/disinfectant (70% alcohol; 500-ppm sodium hypochlorite; hydrogen peroxide, among others). Using a dry paper towel, remove all the disinfectant excess.

The syringe is ready to receive a sterile tip for a new procedure.

⚠ The dental professional and his/her assistant must be wearing IPE during all these steps.

## Maintenance

### Preventive maintenance

In order to reduce equipment faults and increase durability, call Dabi Atlante Authorized Technical Assistance and set a regular preventive maintenance plan.

### Corrective maintenance

In case the equipment presents any faults that are not listed in "Troubleshooting", contact Dabi Atlante Authorized Technical Assistance.

⚠ Do not open the equipment and/or try to repair it as such action may aggravate the problem or even generate further malfunctioning.

### Circuit diagrams, list of parts, components, others

In case they are necessary, ask Dabi Atlante for electric diagrams and/or lists of parts, components or other information, which will be supplied by agreement.

## Flush system, handpieces and jet/ultrasound tank supply (Cuspidor)

Hold the reservoir tightly on its base and give it 1/4 of a turn to the left so as to allow the pressurized air to exit through the grooves. Next, proceed to unscrew it completely and discard the solution remaining on the bottom of the container. Supply the reservoir with new solution and screw it to its support, checking for air leaks.

⚠ The jet/ultrasound solution reservoir is not pressurized.  
Replace the pet tank every 6 months.

## Solutions used

### Flush system

Chlorinated water at 1:400. This concentration can be achieved by dissolving 25 ml of sodium hypochlorite solution at 10,000 ppm in 500 ml of water

### Handpieces

Filtered or distilled water. A weakly chlorinated solution can be used at a concentration of 2 to 3 ppm of chlorine. Dissolve 0.3 ml of sodium hypochlorite at 10,000 ppm in 500 ml of water

### Jet/Ultrasound

Antiseptic solutions such as PVPI (topical), Cepacol, chlorhexidine and Listerine diluted at 50% or more, depending on patient's sensitivity or level of infection to be treated

## Recommendations and care

The use of protection goggles (not supplied with the equipment) is recommended during the operation of Ultralux and Ultraled.

Use dry oil-free filtered air. Install an air filter on you compressor outlet (the filter is not supplied with the equipment). The compressor must be installed at a site where it can be easily handled, since it is necessary to clean every day in case water remains inside it.

Keep the compressed air supply according to the recommended amount and pressure.

Use only the solutions recommended in this manual in the reservoirs.

Perform the cleaning and disinfection processes as recommended.

Make sure that the solution reservoir is not used in the wrong function.

Drain the reservoir solutions at the end of the working day and make the flush system function.

The use of disposable physical barriers is recommended for contact parts such as handles. Transparent PVC film can be used.

## Procedure in case equipment functioning is altered

In case the equipment presents any faults, firstly check whether they are listed under "Troubleshooting". If the problem cannot be solved, turn off the equipment, unplug it from the power source outlet and call for Dabi Atlante Authorized Technical Assistance.

## Sensitivity to predictable environmental conditions in normal working situations

This equipment is not sensitive to magnetic, electric and electrostatic or pressure interferences as long as the cleaning, maintenance, transportation and operation recommendations found in this manual are observed.

## Procedures in case of equipment disposal

When the equipment is put out of use, it must be disposed of at an appropriate site (according to local regulations), thus preventing its inadequate use and environmental contamination.

## Transportation and storage conditions

The equipment must be transported and stored:

- Carefully in order to prevent falls and impact.
- By keeping the side indicated by the arrow facing upwards.
- By observing maximum piling limits as indicated on the packaging.
- Protected against moisture, rain, water sprinkling and wet ground.

## After installation

Keep the equipment protected from sunlight and rain.

Transportation /storage environmental conditions	
Transportation and storage environmental temperature range	0°C to + 55°C
Transportation and storage relative humidity range	0 % to 90 % (non-condensing)
Atmospheric pressure range	500hPa to 1060 hPa (375 mmHg to 795 mmHg)

## Preparation prior to delivery

This equipment has been revised and approved according to the Dabi Atlante Quality Inspection and Control Program for better performance.

Environmental operation conditions	
Environmental temperature functioning range	+10°C to +35°C
Environmental temperature range recommended by Dabi Atlante	+21°C to + 26°C
Functioning relative humidity range	30% to 75 % (non-condensing)
Atmospheric pressure range	700 hPa to 1060 hPa (525 mmHg to 795 mmHg)

## Environmental accommodation operation conditions (between operations)

Environmental temperature accommodation range	+5°C to +45°C
Environmental temperature range recommended by Dabi Atlante	+15°C to +30°C
Functioning relative humidity range	30% to 75 % (non-condensing)
Atmospheric pressure range	700 hPa to 1060 hPa (525 mmHg to 795 mmHg)

## Dabi Atlante Authorized Service Network

Installation and all maintenance service for Dabi Atlante equipment must be given by Dabi Atlante Authorized Technical Assistance, otherwise warranty will be compromised.

Use the attached insert or access the website [www.dabi.com.br](http://www.dabi.com.br) to find the nearest Dabi Atlante Authorized Technical Assistant.

## Symbols

Use the icons below to identify the symbols for you equipment:

### Symbols on the packaging

	<b>"Fragile"</b> Located on the packaging side. It determines that equipment must be carefully transported in order to prevent falls or shocks
	<b>"Protect against moisture"</b> Located on the packaging side. It determines that the equipment must be protected against any type of moisture during transportation and storage
	<b>"This side up"</b> Located on the packaging side. It determines that the packaging must be handled by keeping the side shown by the arrow facing upwards
	<b>"Maximum piling"</b> Located on the packaging side. It determines the maximum number of boxes that can be piled up during transportation and storage
	<b>"Temperature limit"</b> It determines the temperature limits at which the equipment must be stored or transported

### Electronic controls

	Press for the zero position
<b>M1</b>	* Press to program the function
<b>M2</b>	* Press to program the function
<b>M3</b>	* Press to program the function
	Press to activate/turn off the dental light

\* Press to activate/interrupt the water flow in the cuspidor basin

Press to activate/interrupt the water flow in the cup filler

Press to raise the seat

Press to reach the spitting position

Press to lower the seat

Press to lower the backrest

Press to raise the backrest

## Final recommendations

In order to ensure greater durability for your equipment, use only original Dabi Atlante replacement parts, which are in accordance with the technical standards and specifications required by the manufacturer. Also, use only the services provided by the Dabi Atlante Authorized Technical Assistance Network as it counts on trained assistants and specific tools for the appropriate maintenance of your equipment.

## Classification and type

- According to Standard IEC 60601-1: Class I, Type B – Intermittent operation.
- According to Directive 93/42: Class II a.

## Equipment Warranty

For sales outside Brazil, the "Warranty Terms" are the local distributor responsibility.

## Delivery unit technical features (for all models)

Frequency	50/60 Hz
Voltage inside the equipment	24 Vac supplied by the chair
Air pressure	60 to 80 PSI±2
Air consumption	65 l/min
Water consumption	0.1 l/min
Power consumption	140 VA ± 10%
Handpiece connection	Optional 2 or 4 holes according to the ISO 9168 Standard
Load capacity	1.0 Kg
Reservoirs for pressurized solutions	1 500-ml reservoir for handpiece and syringe cooling water – work pressure of 40 psi 1 500-ml reservoir for the flush asepsis system chlorinated water – work pressure of 40 psi

## Curing light technical features

Power consumption	Ultraled: 84 VA ± 10%
Light conductor	11 mm in diameter of semi-oriented fiber optics (code 42842-002/8) or of optional oriented fiber optics (code 26419-032/2)

## Ultrasound technical features

Handpiece vibration frequency	24 to 30 KHz
Handpiece vibration amplitude	From 0.05 to 0.1 mm
Irrigating solution consumption	28 ml/min (with the use of a peristaltic pump)
Power consumed	15 VA ± 10%
Transducer system	Piezo electric ceramic

## Troubleshooting

### Delivery units

Faults	Causes	Solutions
No water in the syringe	No water in the reservoir	Supply reservoir with water
No water or air in the syringe	Compressor air gauge is closed	Open air gauge
No spray in the handpieces	No water in the reservoir	Supply reservoir with water
	Spray adjustment meter is closed	Open meter
Handpieces will not work	Compressor air gauge is closed	Open gauge
Equipment will not go on	Fuse is burned out	Replace fuse
	No power supply	Await normalization
	Circuit breaker is off	Turn circuit breaker on
No pressurization in the system	Compressor is off and/or gauge is closed	Turn compressor on and/or open gauge

### Curing light

Faults	Causes	Solutions
Equipment will not go on	No power supply	Await normalization or call for technical assistance
Equipment will not cure	LED is burned out	Replace head (Ultraled)
	Light conduction tip is blocked	Replace or clean

### Sodium Bicarbonate Jet

Faults	Causes	Solutions
Bicarbonate in the jet is insufficient	Not enough bicarbonate in the reservoir	Supply the reservoir with bicarbonate
	Handpiece tip is clogged	Unclog tip
Pressure in the jet is insufficient	Compressor is off	Turn compressor on
Water in the jet is insufficient	Water gauge is closed	Open gauge
	Water reservoir is empty	Supply reservoir with water

### Ultrasound

Faults	Causes	Solutions
Tip will not oscillate, but light goes on	Tip is not fitted appropriately	Fit tip appropriately or try another tip
Ultrasound handpiece leaks water	Tip is not screwed appropriately	Screw tip appropriately or replace O'ring
Water will not come out	Gauge is closed	Open gauge
Perception of electric shock in the tip	Grounding is inexistent or inefficient	Ground equipment appropriately

**Technical responsibility by:** Leonel Issa Halak – CREA-SP 0600256918

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**MPC International S.A** - 26, Boulevard Royal - L-2449 - Luxembourg - Fone: +49 431-7028 219 - Fax: +49 431-7028 202
MPC International S.A - 26, Boulevard Royal - L-2449