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owner's manual

UNIDADES

cuspidor

DABI ATLANTE

Português.....03
English.....14
Versión en español deste manual:
www.dabi.com.br/esp

Unidade Versátil / Versátil Cuspidor



fig. A

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Contents

Legend	14
Versátil, Versátil Assistant Cuspidor	14
Techno 200, Techno 200 Assistant Cuspidor.....	14
Techno Orto Integrated Cuspidor	15
Installation	15
Procedures prior to equipment use/ reutilization	15
Operation	15
Three-way syringe activation	15
Cuspidor water supply activation	15
Aspirator activation	16
Cup filler activation	16
Ultraled operation	16
Chair activation by the cuspidor panel	17
Light activation by the electronic panel	17
Cleaning	17
Maintenance	19
Recommendations and care	19
Flush system, handpieces and jet/ultrasound tank supply	19
Procedure in case equipment functioning is altered	19
Sensitivity to predictable environmental conditions in normal working situations	19
Procedures in case of equipment disposal	19
Transportation and storage conditions	19
After installation	20
Preparation prior to delivery	20
Dabi Atlante Authorized Service Network	20
Symbols	20
Final recommendations	21
Classification and type	21
Equipment Warranty	21
Cuspidors Technical Features	21
Curing light Technical Features	21
Troubleshooting	22

Presentation

Dear Client,

Welcome!

We are really proud of your preference for our products. Your familiarization with the **Dabi Atlante Dental Cuspidors** will help maximize their performance and increase their durability and benefits. Therefore, before starting your equipment, be sure to read this manual very attentively.

Important

All the information, illustrations and specifications contained in this manual are based on data existing at the time of its publication. We reserve the right to make any modifications deemed necessary both on the product and in this manual without previous notice.

Legend

Versátil Cuspidor (fig. A), **Versátil Assistant Cuspidor** (fig. B)

- ① Basin
- ② Water adjustment tap
- ③ Optional terminal for the Venturi/Vac Plus suction system or high-power (optional)
- ④ Venturi aspirator terminal
- ⑤ Three-way syringe (optional)
- ⑥ Basin-cleaning water pipe
- ⑦ High-impact polystyrene framework
- ⑧ Cup filler
- ⑨ Assisting arm (optional)
- ⑩ Cuspidor control panel (optional)
 - ⑩a Seat ascent
 - ⑩b Backrest ascent
 - ⑩c Light
 - ⑩d Spitting position
 - ⑩e Basin water activation
 - ⑩f Cup filler activation
 - ⑩g Seat descent
 - ⑩h Backrest descent
 - ⑩i Zero position
 - ⑩j Work positions
- ⑪ Electric cuspidor and cup filler water activator
- ⑫ Solution reservoir
- ⑬ Basin drain

Techno 200 (fig. C), **Techno Assistant 200 Cuspidor** (fig. D)

- ① Basin
- ② Water adjustment tap
- ③ Optional terminal for the Venturi/Vac Plus suction system or high-power (optional)
- ④ Venturi aspirator terminal
- ⑤ Three-way syringe (optional)
- ⑥ Basin-cleaning water pipe
- ⑦ High-impact polystyrene framework
- ⑧ Cup filler
- ⑨ Assisting arm (optional)
- ⑩ Cuspidor control panel (optional)
 - ⑩a Seat ascent
 - ⑩b Backrest ascent
 - ⑩c Light
 - ⑩d Spitting position
 - ⑩e Basin water activation
 - ⑩f Cup filler activation
 - ⑩g Seat descent
 - ⑩h Backrest descent
 - ⑩i Zero position

- ⑩ Work positions
- ⑪ Electric cuspidor and cup filler water activator
- ⑫ Solution reservoir
- ⑬ Basin drain

Versátil Orto Integrated Cuspidor (fig. E)

- ① Basin
- ② Water adjustment tap
- ③ Venturi aspirator terminal
- ④ Micromotor terminal
- ⑤ High-speed handpiece terminal
- ⑥ Three-way syringe (optional)
- ⑦ Basin-cleaning water pipe
- ⑧ High-impact polystyrene framework
- ⑨ Assisting arm (optional)
- ⑩ Cuspidor control panel (optional)
 - ⑩a Seat ascent
 - ⑩b Backrest ascent
 - ⑩c Light
 - ⑩d Spitting position
 - ⑩e Basin water activation
 - ⑩f Cup filler activation
 - ⑩g Seat descent
 - ⑩h Backrest descent
 - ⑩i Zero position
 - ⑩j Work positions
- ⑪ Instrument holder
- ⑫ Solution reservoir
- ⑬ Curing Light (optional)
- ⑭ Handle
- ⑮ Cup filler
- ⑯ Basin drain

Installation

This equipment must be unpacked and installed by an authorized Dabi Atlante technician; otherwise, warranty will be voided.

Procedures prior to equipment use/reutilization

Prior to equipment use/reutilization, follow the operation and cleaning instructions in this manual.

Operation

Turning the equipment on/off

Turn the chair's master switch on. All the equipment's functions will be enabled for operation. Turn the master switch off after use.

The master switch features an internal LED that remains on/off when the chair is on/off.

Basin water flow adjustment

Adjust water flow using the adjustment tap (item 2, fig. A, B, C, D and E)

Activation of the Conventional Venturi System aspirators

Remove the aspirator (fig. A, B, C and D item 4 and item 3 fig.E) from its holder and it will automatically start functioning. Place the aspirator back on its holder to turn it off.

High-speed handpiece and micromotor use (Orto Cuspidor)

Connect the hose terminal onto the rear part of the high-speed handpiece or micromotor so that the larger tap hole coincides with the larger hole on the handpiece.


Screw the tap coupling nut onto the handpiece framework (clockwise).

Activate the footswitch and the selected handpiece will start functioning.

The footswitch air outflow is controlled by the operator, and the handpiece rotation can be increased or decreased by applying higher or lower pressure on the footswitch.

Procedure for high-speed handpiece and micromotor spray adjustment (Orto Cuspidor)

Adjustment is made by a meter located on the hose near the terminal. Turn it clockwise to increase spray and counterclockwise to decrease it.

Three-way syringe activation (optional)

Press the green button for water outflow, the blue button for air or both buttons simultaneously for spray.

Cuspidor water supply activation (optional)

Electric activation

Press the switch (fig. A, B, C and D, item 11) for a continuous flow of water into the basin.

Electronic activation

Press key  with a quick touch in order to open or close water flow into the cuspidor.

Press key  with two quick touches to open the water flow into the cuspidor for the programmed time.

Time programming for water flow into the cuspidor's basin

Keep key  pressed and a polyphonic sound will be emitted, indicating the program mode. Release the key after:

- The emission of the first tone to program for 10 seconds;
- The emission of the second tone to program for 20 seconds, successively following the same procedure for up to 60 seconds.

Aspirator activation

Fig.A, B, C, and D, item 4 and fig.E item 3)

Vac Plus and Venturi

Remove the aspirator from its holder and it will automatically start functioning. Place the aspirator back on its holder to turn it off.

High power

Remove the aspirator from its holder and it will be enabled for functioning. Open/close the tap located on the framework of the aspirator under use in order to activate/interrupt suction.

Place the aspirator back on its holder in order to disable it.

Cup filler activation (optional)

Fig. A, B, C and D item 8 and Fig.E item 3

Electric activation

Press the switch (fig. A, B, C and D item 11) to start the water flow.


Electronic activation

Keep key  pressed for the desired time to be programmed. Later on, only a single touch will be necessary to activate it.

Ultraled operation (optional)

Curing

Attach the acrylic curing light tip to the front orifice of the handpiece.

-  The acrylic curing light is supplied with **Ultraled**. It is possible to use the semi-guided tip, which is sold separately. After the tip has been attached to the equipment, it should be cleaned with some gauze, since oiliness from the user's hands and other stains on the acrylic tip may affect light conduction.

Move the handpiece into the patient's mouth and position the light tip at a maximum distance of 2mm from the material to be cured.

Press the button to activate the LED, which will light up slowly. As soon as it is activated, the LED will remain on during the set time and a tone will be heard every 10 seconds.

When the programmed time has expired, the LED will go off.

-  In case curing needs to be repeated, press the Ultraled activation trigger again. In order to interrupt the curing cycle, also press the activation trigger.

Whitening

Attach the acrylic bleaching tip (for 1 or 3 teeth) to the front orifice of the handpiece.

-  The bleaching tips are sold separately. After the tip has been attached to the equipment, it should be cleaned with some gauze, since oiliness from the user's hands and other stains on the acrylic tip may affect light conduction.

Follow the gel manufacturer's instructions as regards activation time and the distance between the tip and the gel.

Press the button to activate the LED, which will light up slowly. As soon as it is activated, the LED will remain on during the set time and a tone will be heard every 10 seconds.

When the programmed time has expired, the LED will go off. (Refer to Time setting).

-  To interrupt the bleaching cycle, press the activation button again.

Time Setting

Ultraled

Press the activation button on the handpiece to activate the LED. Then, without releasing the activation button, wait three seconds, and a long tone will be emitted. Still keeping the activation button pressed, wait until a second tone is emitted. Release the activation button and wait for a tone indicating the first setting phase, "Low Power". Then, press the activation button to set the application percentage. Each touch corresponds to 10% of the total preadjusted total functioning time, and the maximum will be 100% (in case the activation button is not pressed, the system will understand that no percentage has been set in this phase).

After the desired application percentage has been programmed, and if it is less than 100%, wait 2 seconds, and 2 tones will be emitted, indicating the second setting phase, "Medium Power". Repeat the procedure used for the first phase. At this moment, only the remaining percentage can be programmed. "High Power" will correspond to the remaining percentage from the two first phases. Two seconds after the setting of the second phase has been concluded, 3 fast tones will be emitted, indicating setting completion. In case you wish to set only some of the phases, access the setting model as explained above, but do not press the activation button after each phase. The system will automatically shift to the next phase. When the activation button is not pressed during the setting of phase 1 and 2, the application will be set at 100% with "High Power". Example: For an application time setting of 60 seconds and the phases respectively set at 10, 20

and 70%, the LED will remain on for 6 seconds at low power, for 12 seconds at medium power and for 42 seconds at high power.

Chair activation by the cuspidor panel (optional)

Sound indication configuration

1. Place the chair at the zero position, by pressing the zero-position key (fig. A, B, C, D and E, item 10i);

2. Keep the zero-position key pressed. A long sound signal will be emitted, indicating the setting mode. Release the key after:

- The first tone (single) for configuration without a sound signal;
- The second tone (double) for configuration with a standard sound signal;
- The third tone (triple) for configuration with a polyphonic sound signal.

 Regardless of the configuration selected, the work position controls will be followed by a sound signal, and the movements activated by the joystick will only emit a sound signal when the seat and backrest reach their highest or lowest position.

Seat movement

Press key  /  to raise or lower the seat.

Backrest movement

Press key  /  to raise or lower the backrest.

Automatic movement to zero position

Press key  for the seat and backrest to simultaneously return to the zero position. This is the most comfortable position for the patient to get on/off the chair.

Automatic movement to the work positions

Press and release key:

- **1** to move the chair to work position 1;
- **2** to move the chair to work position 2;
- **3** to move the chair to work position 3.

After releasing the key, the chair's seat and backrest will simultaneously move to the selected work position and the light will be set to the programmed intensity.

In order to interrupt the movement, activate any other control.

Work positions programming

1. Place the chair at the desired work position using the joystick;
2. Keep keys 1, 2 or 3 pressed to save the chair's position in memory **1**, **2** or **3**, respectively.

 It is possible to save the light's intensity and the work position at the same time, by setting it to the desired intensity prior to beginning the saving process.

Automatic movement to the spitting position

In order to move the backrest to the spitting position, quickly press key .

In order to resume the previous position, quickly touch key  twice.

In order to interrupt the movement, quickly touch the spitting position key while the movement is being performed.

Light activation by the electronic panel (optional)

Press and release key  to turn the light on or off.

 The reflector must be connected to the chair's circuit.

In order to change light intensity with the **Croma** and **Galla** chairs with a Fuse system, keep the key pressed until the desired intensity is achieved.

Flush system, handpieces and jet/ultrasound tank supply

Hold the reservoir tightly on its base and give it $\frac{1}{4}$ of a turn to the left so as to allow the pressurized air to exit through the grooves. Next, proceed to unscrew it completely and discard the solution remaining on the bottom of the container. Supply the reservoir with new solution and screw it to its support, checking for air leaks.

 The jet/ultrasound solution reservoir is not pressurized. Replace the PET reservoirs every six months.


Cleaning

Clean the equipment surface with a clean soft cloth moistened in Aplic Odonto® by Dabi Atlante or another product with similar chemical properties.

 The use of other chemical products is not recommended as it may damage the equipment.

Cleaning the cuspidor's basin


Remove the Techno Cuspidor's basin.


Wash it using a soft sponge, water and neutral soap. Before returning the basin to its place, lubricate the rubber ring with Vaseline in order to prevent it from drying up, which will make removal and placement easier.

Cleaning the waste separator

Remove the lid of the waste separator located on the chair's cover by turning it counterclockwise.

The waste separator's screen will be removed together with the lid. It must then be washed in running water.

Return it to place by turning it clockwise.

Cleaning aspirator ducts

Venturi aspirator

Prepare a solution by adding 30ml of germicide to 500ml of water and aspirate 250ml of such solution with both aspirators.

Pour the remaining 250ml into the basin.

Aspirate 250ml of clean water with both aspirators.

Let water flow into the basin for 1 to 2 minutes.

This procedure must be performed on a daily basis after the working hours.

High-power/Vac Plus aspirator

Prepare a solution by adding 60ml of germicide to 1000ml of water and aspirate 750ml of such solution with the aspirators.

Pour the remaining 250ml into the basin.

Aspirate 500ml of clean water with the aspirators.

Let water flow into the basin for 1 to 2 minutes.

This procedure must be performed on a daily basis after the working hours.

Sterilization of the high-power suction electors


Wash the handpieces in running water and using neutral soap in order to remove external and internal residues.

Rinse and dry the handpieces thoroughly prior to placing them in the autoclave.


The blood and saliva ejectors are autoclavable.

Lubrication of the high-power suction ejectors.

After cleaning and sterilizing the ejectors, remove the activation lever by supporting the ejector on a piece of furniture and firmly pressing the lever with your thumb.


Remove the ejector's cylinder and wash it in running water.

Lubricate the cylinder's rings with Vaseline and replace it on the elector.

Replace the activation lever on the ejector by pressing firmly.

Replace the ejector onto the cuspidor's hose connection.

Failure in cleaning and lubricating will cause permanent damage to the ejector.

Cleaning the curing light

Clean the surface of the Ultraled or Ultralux with a clean soft cloth moistened in Aplic Odonto® by Dabi Atlante or another product with similar chemical properties. For better preservation of tips, prevent accidents that may cause damage and try to store them in a dry and airy environment.


In case PVC film is used, make sure that the ventilation orifices (air inlet and outlet) have not been obstructed. Never wrap the acrylic tip with any type of material.

Sterilization of the curing light tips


The acrylic curing light tip as well as the bleaching light tips must be sterilized by cold sterilizing solution and must not be autoclaved.

Cleaning the terminals support cover

Remove the handpieces from their supports.
Hold the cover and pull the locks with your finger in order to release them.

Wash the cover using a soft sponge, water and neutral soap
Replace the cover following the inverse procedure.

Make sure that the cover is securely placed.
Return the handpieces to the supports.

Maintenance

Preventive maintenance

In order to reduce equipment faults and increase durability, call Dabi Atlante Authorized Technical Assistance and set a regular preventive maintenance plan.

Corrective maintenance

In case the equipment presents any faults that are not listed in "Troubleshooting", contact Dabi Atlante Authorized Technical Assistance.

Do not open the equipment and/or try to repair it as such action may aggravate the problem or even generate further malfunctioning.

Circuit diagrams, list of parts, components, others

In case they are necessary, ask Dabi Atlante for electric diagrams and/or lists of parts, components or other information, which will be supplied by agreement.

Recommendations and care

The use of protection goggles (not supplied with the equipment) is recommended during the operation of Ultraled.

Use dry oil-free filtered compressed air. Install an air filter on you compressor outlet (the filter is not supplied with the equipment). The compressor must be installed at a site where it can be easily handled, since it is necessary to clean every day in case water remains inside it.

Keep the compressed air supply according to the recommended amount and pressure.

Use only the solutions recommended in this manual in the reservoirs.

Perform the cleaning and disinfection processes as recommended.

Drain the reservoir solutions at the end of the working day and make the flush system function.

The use of disposable physical barriers is recommended for contact parts such as handles. Transparent PVC film can be used.

Procedure in case equipment functioning is altered

In case the equipment presents any faults, firstly check whether they are listed under "Troubleshooting". If the problem cannot be solved, turn off the equipment, unplug it from the power source outlet and call for Dabi Atlante Authorized Technical Assistance.

Sensitivity to predictable environmental conditions in normal working situations

This equipment is not sensitive to magnetic, electric, electrostatic and pressure interferences as long as the cleaning, maintenance, transportation and operation recommendations found in this manual are observed.

Procedures in case of equipment disposal

When the equipment is put out of use, it must be disposed of at an appropriate site (according to local regulations), thus preventing its inadequate use and environmental contamination.

Transportation and storage conditions

The equipment must be transported and stored:

- Carefully in order to prevent falls and impact.
- By keeping the side indicated by the arrow facing upwards.
- By observing maximum piling limits as indicated on the packaging.
- Protected against moisture, rain, water sprinkling and wet ground.

Transportation /storage environmental conditions	
Transportation and storage environmental temperature range	0°C to + 55°C
Transportation and storage relative humidity range	0% to 90% (non-condensing)
Atmospheric pressure range	500hPa to 1060 hPa (375 mmHg to 795 mmHg)

After installation

Keep the equipment protected from sunlight and rain.

Environmental operation conditions	
Environmental temperature functioning range	+10°C to +35°C
Environmental temperature range recommended by Dabi Atlante	+21°C to + 26°C
Functioning relative humidity range	30% to 75% (non-condensing)
Atmospheric pressure range	700 hPa to 1060 hPa (525 mmHg to 795 mmHg)

Environmental accommodation operation conditions (between operations)	
Environmental temperature accommodation range	+5°C to +45°C
Environmental temperature range recommended by Dabi Atlante	+15°C to +30°C
Functioning relative humidity range	30% to 75% (non-condensing)
Atmospheric pressure range	700 hPa to 1060 hPa (525 mmHg to 795 mmHg)

Preparation prior to delivery

This equipment has been revised and approved according to the Dabi Atlante Quality Inspection and Control Program for better performance.


Dabi Atlante Authorized Service Network

Installation and all maintenance service for Dabi Atlante equipment must be given by Dabi Atlante Authorized Technical Assistance, otherwise warranty will be compromised. Use the attached insert or access the website www.dabi.com.br to find the nearest Dabi Atlante Authorized Technical Assistant.

Symbols

Use the icons below to identify the symbols for you equipment:

Symbols on the packaging

	"Fragile" – Located on the packaging side. It determines that equipment must be carefully transported in order to prevent falls or shocks.
	"Protect against humidity" – Located on the packaging side. It determines that the equipment must be protected against any type of humidity during transportation and storage.
	"This side up" – Located on the packaging side. It determines that the packaging must be handled by keeping the side shown by the arrow facing upwards.
	"Maximum piling" – Located on the packaging side. It determines the maximum number of boxes that can be piled up during transportation and storage.
	"Temperature limit" – It determines the temperature limits at which the equipment must be stored or transported.

Electronic controls

	Press for zero position
	Press to program work position
	Press to program work position
	Press to program work position
	Press to activate/turn off the dental light
	Press to activate/interrupt the water flow in the cuspidor basin
	Press to activate/interrupt the water flow in the cup filler
	Press to raise the seat
	Press to lower the seat
	Press to lower the backrest
	Press to raise the backrest
	Press for the spitting position

Final recommendations

In order to ensure greater durability for your equipment, use only original Dabi Atlante replacement parts, which are in accordance with the technical standards and specifications required by the manufacturer. Also, use only the services provided by the Dabi Atlante Authorized Technical Assistance Network as it counts on trained assistants and specific tools for the appropriate maintenance of your equipment.

Classification and type

- According to Standard IEC 60601-1: Class I, Type B – Intermittent operation.
- According to Directive 93/42: Class II a.

Equipment Warranty

The “Warranty Terms” are found on the certificate accompanying the product. Any alterations or concessions involving the warranty as well as pertaining conditions and authorizations, either oral or written, are definitely not permitted if not previously documented and approved of by Dabi Atlante.

Cuspidors Technical Features

Mains	Supplied by the chair
Voltage	24 Vac
Frequency	50/60 Hz
Power consumption	-
Timer and cup holder activation	20VA±10%

consumption	Air consumption	W a t e r l/min l/min
Vac Plus	35	-
Venturi	9	-
Cuspidor with timer/curing light	35	2
Micromotor	65	0.02
High-speed handpiece	40	0.02
Three-way syringe	17	0.1

Curing light Technical Features

Especifications	Ultraled
Voltage	Comutação Automática 100 a 240 VCA
Consumption	10 VA
Frequency	50/60 Hz
Length	205 mm
Width	138 mm
Height	200 mm
Weight	0,65 Kg

Troubleshooting

Faults	Causes	Solutions
Weak suction by the aspirator	Waste accumulation in the separator's screen on the hose line	Clean the waste separator's screen Open the compressor's air tap
Aspirator will not work	No air supply in the system	Await normalization
No water flow in the cuspidor	No water supply or cuspidor tap is closed	Open the water tap Await normalization
Curing light will not go on	No power supply	Replace head
Equipment will not cure	LED is burned out	Replace tip or clean it
	Light conduction tip is clogged	Replace bulb
	Bulb is burned out	

Technical responsibility by: Leonel Issa Halak – CREA-SP 0600256918

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MPC International S.A - 26, Boulevard Royal - L-2449 - Luxembourg - Phone: +49 431-7028 219 - Fax: +49 431-7028 202
MPC International S.A - 26, Boulevard Royal - L-2449